

MADE IN HUNGARY • PI DENTAL MANUFACTURING CO. LTD.

# CYBER SCAN<sup>®</sup>

3D Dental Scanner

*PI dental*<sup>®</sup>  
PIDENTAL.HU



High quality  
at reasonable price

Compatible with  
**exocad**

# CYBER-SCAN®

## CAD/CAM solutions for dental labs and dentists



### TECHNICAL DATA


Main supply	230 V, 50/60 Hz
Axis	2 axis
Light Source	White light LED
Technology	Structured light
Camera Resolution	1,3 MP
Software	Open Type (STL compatibility)
Accuracy	10 µm
Weight	16 kg
O/S	Windows 7/8/10 (64bit)
Dimension (Width x Height x Depth)	310 x 440 x 330 mm

- Optimized for productive, fast and accurate scanning
- High productivity
- 2 cameras with 1,3 megapixels
- Fully compatible with EXOCAD®
- Order No.: 77-0001

# EXOCAD®


### FOR DENTAL TECHNICIANS

- Be part of the future. EXOCAD® DentalCAD is the modern tool for every dental technician who wants to extract the advantages from the ongoing digital dental laboratory revolution.
- To digitally design a restoration has never been as easy and accessible - without the steep learning curve of the past. This will help you maximize productivity and profits. EXOCAD® DentalCAD is reliable and robust even when dealing with complex cases on a daily basis.

**exocad**

### HIGH TECH TOOLS FOR HIGH PRECISION ARTISANS

- Already the standard version our software covers a wide variety of indications, making it an attractive choice both for expert users and newcomers in the dental CAD field. The strenght of EXOCAD® DentalCAD really lies in its simplicity. There is no point in creating software that makes your job harder, so it doesn't. Even for first time users, it's an extremely logical and easy to use program that emulates the way technicians work in the analog world.
- This software grows with your needs: for advanced indications, such as implantology, a variety of add-on modules is available.



• Order No.: 78-0000-01

**2017**

Published with the support of:



HUNGARIAN  
INVESTMENT AND  
TRADE AGENCY

**PI DENTAL FOGÁSZATI GYÁRTÓ KFT. (MANUFACTURING CO. LTD.)**

H-1141 Budapest, Szugló u. 83-85., Hungary • Phone: +36 1 2514944, +36 1 3632234  
Fax: +36 1 2514891 • E-mail: sales.pidental@pidental.hu • www.pidental.hu