

MADE IN HUNGARY • PI DENTAL MANUFACTURING CO. LTD.

CYBER SCAN ART PLUS®

3D Dental Scanner



High quality
at reasonable price

Compatible with
exocad

CYBER SCAN ART PLUS®

CAD/CAM solutions for dental labs and dentists



- Time saving scan process
- No modell fixture
- Direct articulator scanning
- Multi-die Scanning
- Implant Multi-die Scanning
- Impression Scanning
- Orthodontic Application
- Pencil Mark Scanning
- All in one scan
- Fast scanning: 27 sec SINGLE DIE
50 sec FULL MODEL
- Fully compatible with EXOCAD®
- ORDER No.: 77-0002

TECHNICAL DATA

Main Supply	230 V, 50/60 Hz
Axis number	2
Light Source	White light LED
Technology	Structured light
Camera Resolution	2,0 MP
Software	Open Type (EXOCAD, 3shape, STL compatibility)
Accuracy	10 µm
Dimension (Width x Height x Depth)	330 x 495 x 430 mm
Weight	15 kg

Scanning speed

Speed (one die scan)	27 sec
Speed (Full arch)	50 sec

Scanning methods

Impression scanning	Yes
Multi-die scanning	Yes
Texture scanning	Yes
Implant multi-die scanning	Yes
Articulator scanning	Yes (type A,B,S)

System requirement

Operating system	Windows 7/8/10 (64-Bit)
RAM	16 GB DDR3
Video Card	NVIDIA GT750Ti (2GB)
HDD	500 GB
CPU	Intel Core i7 or i5
Chipset	Z97, H97, X99
USB port	USB 3.0


EXOCAD®

FOR DENTAL TECHNICIANS

- Be part of the future. EXOCAD® DentalCAD is the modern tool for every dental technician who wants to extract the advantages from the ongoing digital dental laboratory revolution.
- To digitally design a restoration has never been as easy and accessible - without the steep learning curve of the past. This will help you maximize productivity and profits. EXOCAD® DentalCAD is reliable and robust even when dealing with complex cases on a daily basis.



exocad



• Order No.: 78-0000-01

2017

Published with the support of:



HUNGARIAN
INVESTMENT AND
TRADE AGENCY



PI DENTAL FOGÁSZATI GYÁRTÓ KFT. (MANUFACTURING CO. LTD.)

H-1141 Budapest, Szugló u. 83-85., Hungary • Phone: +36 1 2514944, +36 1 3632234
Fax: +36 1 2514891 • E-mail: sales.pidental@pidental.hu • www.pidental.hu