


ANGULATED DAS MULTI-UNIT



DIFFERENT **GINGIVAL HEIGHTS**


20° ANGULATION

4MM PLATFORM DIAMETER

SAME COMPONENTS AS THE STRAIGHT
MULTI-UNIT

DYNAMIC SCREWDRIVER IS USED
FOR ITS PLACEMENT

ANGULATED DAS MU



ANGULATED DAS MU INSERTION TOOL

COMPATIBILITY	GH (mm)			
	1,5	2,5	3,5	4,5
0004		■	■	■
0006		■	■	■
0007	■	■	■	
0014			■	■
0015	■	■		■
0021	■	■		■
0022	■	■		
0030	■	■	■	■
0033		■	■	■

COMPATIBILITY	GH (mm)			
	1,5	2,5	3,5	4,5
0035		■	■	
0040		■	■	■
0102		■		
0119			■	
0120-0121		■	■	
0125	■	■	■	
0169	■	■	■	
0186		■	■	
0207		■	■	■

CONTACT US FOR MORE COMPATIBILITIES